

 Department of Microbiology Quality Manual	Policy # MI_WL_MSS	Page 1 of 18
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	Subject Title: <b>WASPLab Manual Streaking Station Operator Manual</b>	
Prepared by QA Committee		
Issued by: Laboratory Manager	Revision Date: 1/16/2024	
Approved by Laboratory Director: Microbiologist-in-Chief	Next Review Date: 1/16/2026	

**Uncontrolled When Printed**

### Table of Contents

<b>Introduction.....</b>	<b>2</b>
<b>Manual Streaking Station Installation.....</b>	<b>2</b>
<b>Manual Streaking Station Procedure.....</b>	<b>3</b>
<b>Printer Maintenance .....</b>	<b>6</b>
<b>How to Change the Label Roll.....</b>	<b>6</b>
<b>How to Change the Ink Ribbon .....</b>	<b>12</b>
<b>Reference .....</b>	<b>17</b>
<b>Record of Edited Revisions .....</b>	<b>18</b>

  Department of Microbiology Quality Manual	Policy # MI_WL_MSS	Page 2 of 18
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: WASPLab Manual Streaking Station Operator Manual</b>	

## Introduction

The Manual Streaking Station (MSS) is a WASPLab component used for the manual seeding of Media Plates to be incubated in WASPLab system.

The MSS software communicates with the LIS to provide the operator instructions about streaking of the plates. Those plates are labelled using the MSS printer and are automatically checked in WASPLab for incubation conditions.

This section describes the normal operation of the WASPLab Manual Streaking Station.



## Manual Streaking Station Installation

Before use of the Manual Streaking Station it is necessary to setup the streaking protocols and details in the WASPcore system, this operation is reserved to authorized personnel.

The MSS can work only with an active LIS connection.

The PC where the Manual Streaking Station software is installed must be connected to a printer, to create labels for the streaked media plate.

For details refer to the vendor manual.

  Department of Microbiology Quality Manual	Policy # MI_WL_MSS	Page 3 of 18
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: WASPLab Manual Streaking Station Operator Manual</b>	

## Manual Streaking Station Procedure

1. On the workstation PC desktop, if not already open, double click on the icon “**Manual Streaking Station**”. It may takes approximately 10 seconds for the software to start up.
2. Select the field **Sample barcode**, scan the barcode of the sample and click on **Query**.





UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.*

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_WL_MSS	Page 4 of 18
	Version: 1.1 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>WASPLab Manual Streaking Station Operator Manual</b>	

3. The software screen will population with information about the plates to streak:



- |  |   |
|--|---|
| A. <b>Protocol:</b> protocol name                          | G. <b>Print:</b> print label for one plate  |
| B. <b>Media:</b> progressive number of the plate to streak | H. <b>Print all:</b> print labels for all plates                                  |
| C. <b>Media:</b> type of plate                             | I. <b>First :</b> go back to the first media of the protocol.                     |
| D. <b>Loop:</b> type of loop                               | J. <b>Previous:</b> go back to the previous media                                 |
| E. <b>Streak:</b> the streaking pattern                    | K. <b>Next:</b> go to the next media  |
| F. <b>Media Barcode:</b> barcode of the streaked media     | L. <b>Send to LIS:</b> to send to the LIS the confirmation of the media streaking |



- Click on Print All (H) to print the labels to be applied to the plates
- Streak the plate according to the protocol (A)
- If streaking one plate at a time, click on Next (9), to display the streaking details of the next plate.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed.***

*Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*


Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

  Department of Microbiology Quality Manual	Policy # MI_WL_MSS	Page 5 of 18
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: WASPLab Manual Streaking Station Operator Manual</b>	

- When all the plates for the sample are streaked and labelled, click on **Send to LIS (L)** to send to the LIS the confirmation of the order.

**ATTENTION:** The ‘**Send to LIS**’ button become active only if all the media label are printed.

This command allows the check-in of the plates in WASPLab, that can be manually loaded to start the incubation accordingly to the defined protocol. See [WASPLAB Operator Manual](#).


 Department of Microbiology Quality Manual	Policy # MI_WL_MSS	Page 6 of 18
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: WASPLab Manual Streaking Station Operator Manual</b>	

### *Printer Maintenance*


The Manual Streaking station only requires printer maintenance including replacement of labels and ribbon upon use of the rolls.

### **How to Change the Label Roll**

Below is reported the scheme to be follow for the label roll replacement.




1. Turn lever counterclockwise to unlock the Printhead.




UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.  
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

 Department of Microbiology Quality Manual	Policy # MI_WL_MSS	Page 7 of 18
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title:</b> WASPLab Manual Streaking Station <b>Operator Manual</b>	

- Turn the ring counterclockwise, so that the arrows points to the unlocked “” symbol and releases the margin stop from the roll retainer.



- Pull the margin stop out of the roll retainer.




UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY


*NOTE: This document is **Uncontrolled When Printed**.  
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\




 Department of Microbiology Quality Manual	Policy # MI_WL_MSS	Page 8 of 18
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: WASPLab Manual Streaking Station Operator Manual</b>	

- Remove the previous label roll but leave the black hub on the knob.



- Load the new label roll on the black hub on the knob, and position it below the roller and guide it through the print unit.




UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.  
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\



 Department of Microbiology Quality Manual	Policy # MI_WL_MSS	Page 9 of 18
	Version: 1.1 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>WASPLab Manual Streaking Station Operator Manual</b>	

6. Guide material through the print unit below the label sensor.




7. Move guides closely to the edges of the material without clamping the material.



UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.  
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

 <p>Department of Microbiology</p> <p>Quality Manual</p>	<p>Policy # MI_WL_MSS</p> <p>Version: 1.1 CURRENT</p>	<p>Page 10 of 18</p>
<p><b>Section:</b> Bacteriology Procedures</p>	<p><b>Subject Title:</b> WASPLab Manual Streaking Station Operator Manual</p>	

8. Re-mount the margin stop and push against the label roll as far as possible.




9. Turn ring clockwise, so that the arrow points to the “🔒” symbol, and thus fix the margin stop on the roll retainer.



UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.  
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

 Department of Microbiology Quality Manual	Policy # MI_WL_MSS	Page 11 of 18
	Version: 1.1 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>WASPLab Manual Streaking Station Operator Manual</b>	

10. Make sure that one or two labels are outside the Printhead. This will avoid problems of stuck.




11. Press print head retainer down and turn lever clockwise to lock the printhead.



UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY


*NOTE: This document is **Uncontrolled When Printed**.  
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\


 Department of Microbiology Quality Manual	Policy # MI_WL_MSS	Page 12 of 18
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: WASPLab Manual Streaking Station Operator Manual</b>	

## How to Change the Ink Ribbon

Below is reported the scheme to be follow for the ribbon roll replacement.



1. Turn lever counterclockwise to open the printhead.




UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.  
 Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\




 Department of Microbiology Quality Manual	Policy # MI_WL_MSS	Page 13 of 18
	Version: 1.1 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>WASPLab Manual Streaking Station Operator Manual</b>	

- Remove the used set for transfer ribbon from the two hubs.




- Get the set for transfer ribbon replacement, which is composed by a transfer roll and a ribbon core.

 Department of Microbiology Quality Manual	Policy # MI_WL_MSS	Page 14 of 18
	Version: 1.1 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>WASPLab Manual Streaking Station Operator Manual</b>	



4. Load the transfer ribbon on the hub and the ribbon core on the take-up hub over the printhead.




5. Guide transfer ribbon through the printing unit.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

---

*NOTE: This document is **Uncontrolled When Printed**.*  
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.  
Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

 Department of Microbiology Quality Manual	Policy # MI_WL_MSS	Page 15 of 18
	Version: 1.1 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>WASPLab Manual Streaking Station Operator Manual</b>	



6. Guide transfer ribbon over the label sensor!




7. Secure starting end of the transfer ribbon to the ribbon core using adhesive tape.

UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.  
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\



 Department of Microbiology Quality Manual	Policy # MI_WL_MSS	Page 16 of 18
	Version: 1.1 CURRENT	
<b>Section:</b> Bacteriology Procedures	<b>Subject Title: WASPLab Manual Streaking Station Operator Manual</b>	




- Turn transfer ribbon take-up hub counterclockwise to smooth out the feed path of the transfer ribbon.



UNIVERSITY HEALTH NETWORK/MOUNT SINAI HOSPITAL, DEPARTMENT OF MICROBIOLOGY

*NOTE: This document is **Uncontrolled When Printed**.  
Any documents appearing in paper form that do not state "CONTROLLED COPY" in red print are not controlled and should be checked against the document (titled as above) on the server prior to use.*

Management System\UHN\_Mount Sinai Hospital Microbiology\Standard Operating Procedures\Bacteriology Procedures\

 Department of Microbiology Quality Manual	Policy # MI_WL_MSS	Page 17 of 18
	Version: 1.1 CURRENT	
Section: Bacteriology Procedures	Subject Title: <b>WASPLab Manual Streaking Station Operator Manual</b>	

9. Press printhead down and turn lever clockwise to lock printhead.



### Reference

Manual Streaking Station Operator Vendor Manual, Document number: HPAMSSSEN REV00, COPAN WASP S.r.l.

